

3. Applications

The nomag IC 1.5 MR-Incubator is intended for use for the following purpose in connection with magnetic resonance application; the test certificate refers to this purpose:

1. MR Diagnostics Incubator

4. Restrictions

During the use of this product the following restrictions of functionality and/or applications of the MRI systems, options and accessories apply:

1. Silicone tubes are visible in MR applications with very short TE echo times. A corresponding warning note in the operation instruction is required. Replacement of the Silicone tubes is strongly recommended.

5. Warnings

The following safety precautions must be observed when using this product in connection with the MRI systems listed in section 2:


1. Silicone tubes are visible in MR applications with very short TE echo times. A corresponding warning note in the operation instruction is required.
2. The voltage transformer of the LMT nomag IC 1.5 MR Incubator must be operated at a safe distance of 2m from the magnet and be equipped with a warning sign.
3. The LMT nomag TR must be placed outside the RF-cabin while scanning is taking place.
4. No steel cylinders may be transported with the LMT nomag TR. It and its instruction manual must bear a warning sign.

6. Validity:

This test certificate shall be valid until revoked by Siemens AG, Medical Solutions.

Erlangen, 2007-06-04

Siemens AG



Dr. Söldner



Schumann